

# Lab-Tips<sup>®</sup> - LTP125 Large Knitted Polyester Swab

Cleanroom laundered double layer 100% polyester knit swab

The LTP125 large, knitted polyester swab is an excellent choice for general purpose cleaning or surface sampling. Its durable flat head and rigid internal paddle provide maximum surface area for cleaning or sampling broad flat surfaces. Its cleanroom laundered knitted head provides ultra-low NVR and ion levels.



## **Key Attributes**

- Cleanroom laundered 100% continuous filament polyester knit head
- Long polypropylene handle with a rigid paddle head
- Thermally bonded head
- Laundered and packaged in Class 10 {ISO Class 4} cleanroom

#### Benefits

- · Critically low particles, fibres, ions and extractables
- High abrasion resistance while soft and non-abrasive
- Excellent chemical resistance
- No contaminating adhesives and good sorbency
- Autoclavable
- Free from Silicone, Amides and DOP

### Applications

- Designed for the highest level of contamination control in critical processing applications
- Cleaning of sensitive equipment such as print heads
- Sampling surfaces for cleaning validation
- · Applying lubricants and other liquids
- Removing excess materials
- Cleaning with solvents

## Other Knitted Polyester tip swabs

- LTP1465 Long handle
- LTP 70F Short handle flexible head
- LTP70R Short handle rigid head

## **Physical Characteristics**

| Head Material   | Knitted polyester | Handle Thickness  | 2.8mm {0.110"} |  |
|-----------------|-------------------|-------------------|----------------|--|
| Head Width      | 12.7mm {0.50"}    | Handle Length     | 100mm {3.94"}  |  |
| Head Thickness  | 5mm {0.196"}      | Total Swab Length | 125mm {4.92"}  |  |
| Head Length     | 25mm {0.984"}     | Head Bond         | Thermal        |  |
| Handle Material | Polypropylene     | Handle Text       | Berkshire      |  |
| Handle Width    | 5.5mm {0.216"}    |                   |                |  |

## www.berkshire.uk.com

Contact: Tel + 44 1953 562800 enquiries@berkshire.uk.com

| America | Tel 1 413 528 2602   | info@berkshire.com         |
|---------|----------------------|----------------------------|
| Europe  | Tel + 44 1953 562800 | enquiries@berkshire.uk.com |
| SE Asia | Tel 65 6252 4313     | enquiries@berkshire.com.sg |
| Japan   | Tel 81 3 4530 9883   | master@berkshire.co.jp     |



### **Technical Data:**

| ATTRIBUTE            |             | LTP125 | TEST METHOD               |
|----------------------|-------------|--------|---------------------------|
| Non-Volatile Residue | IPA mg/swab | 0.075  | IEST-RP-CC004.4 Sec 8.1.2 |
| lons (µg/swab)       | Chloride    | 0.11   | IEST-RP-CC004.4 Sec 8.2.2 |
|                      | Sulfate     | 0.061  |                           |
|                      | Nitrate     | ND*    |                           |
|                      | Phosphate   | ND*    |                           |
|                      | Fluoride    | ND*    |                           |
|                      | Potassium   | 0.010  |                           |
|                      | Calcium     | 0.082  |                           |
|                      | Sodium      | 0.032  |                           |
|                      | Magnesium   | 0.011  |                           |

#### Notes:

- Data shown are typical values
- Due to differences in test methods applied and equipment utilised by different swab manufacturers, valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions.
- \* ND\* = below the detection limit of the equipment.

#### **Order Information:**

| Product                 | Item Number | Size  | Swabs/pk | Pks/cs | Case size         |
|-------------------------|-------------|-------|----------|--------|-------------------|
| Large Knitted Polyester | LTP125.5    | 125mm | 100      | 5      | 220 x 160 X 105mm |



## www.berkshire.uk.com

Contact: Tel + 44 1953 562800 / info@berkshire.uk.com